

Dry Mount Tissue

Technical Data

For use in all types of dry mounting presses. For bonding RC photographic paper and printed materials to smooth surfaces.

Arista Dry Mount Tissue is a 50gsm, pH neutral, off-white mounting tissue. It is coated on both sides with a heat activated solvent acrylic adhesive. Available in various sizes.

SPECIFICATIONS:

• Adhesive: Heat Activated Solvent Acrylic

Weight: 50gsmThickness: 3 mil

MATERIAL REQUIREMENTS (for use in a Dry Mount Press):

• Temperature: 185°F to 212°F (85°C to 100°C)

• Dwell Time: 30 to 90 seconds

DISCLAIMER: This is a PERMANENT BOND tissue and all recommendations and instructions are intended as a guide. Users should determine, prior to use, the suitability of each material for their specific purpose. Arista is not responsible for improper use of the product, either by application or substrate to which it is applied.